

## Advanced Technology for Precise Treatments

### Equipped for precision<sup>1</sup>

- Range of exceptional treatment options, including **Contoura<sup>®</sup> Vision** and **StreamLight<sup>™</sup> Transepithelial PRK**
- Enhanced repetition rate of 500 Hz
- Advanced eye-tracking technology
- Non-contact online pachymetry
- High-definition optics
- Integrated **WaveLight<sup>®</sup> Slit Illumination System**



WaveLight® EX500 Excimer Laser

18.08.1973 ♂

OD  Wavefront Optimized

Patient F5

Planning F7

**Treatment F8**

Documentation F9

Setup F10

Laser F11

### REFRACTIVE AND CORNEAL DETAILS

Refraction **-0.50 D -0.25 D x 10°** 12.0 mm

Pupillometry **6.5 mm**

Pachymetry **s 525 t 570 c 550 n 530 i 525 μm**

K1 / Q1 **44.29 D @ 173° / -0.26**

K2 / Q2 **45.06 D @ 83° / -0.26**

### TREATMENT DETAILS

Target Refraction **+0.00 D +0.00 D x 10°** 12.0 mm

Correction **-0.50 D -0.25 D x 10°** 12.0 mm

Optical Zone **6.50 mm** Flap Thickness **120 μm**

Ablation Zone **9.00 mm** Cornea **525 μm**

Residual Stroma **393 μm**

Max. Ablation **12 μm** Centr. Ablation **12 μm**

### ABLATION PROFILE



max **11.67 μm**

**i** Your comment

### TREATMENTS

Patient	Eye	Method	Date	Planned by
Mustermann, Peter	OD	WFO	14.08.2017 11:47:23	Surgeon, Sur

Planned Abort

User friendly graphical User Interface for an enhanced procedural experience



Illuminated keyboard for optimized visibility



Integrated emergency stop for controlled activation

# COLORS TO COMPLEMENT YOU



Custom Colors

---

## Footnotes:

## References

\* As compared to conventional PRK.

+ The results presented in this section are from a prospective, non-randomized study of 230 eyes that had myopia with or without astigmatism and that were treated with topography-guided LASIK with the ALLEGRETTO WAVE<sup>®</sup> Eye-Q Excimer Laser System. Eyes had nearsightedness up to -9 D and astigmatism up to 6 D.

? Pentacam is a trademark of OCULUS Optikgeräte GmbH.

1. Date on File / WaveLight<sup>®</sup> EX500 Excimer Laser User Manual.
2. Mrochen M, Donitzky C, Wüllner C, Löffler J. Wavefront Optimized<sup>®</sup> ablation profiles: theoretical background. J Cataract Refract Surg. 2004;30:775-785.

3. Kanellopoulos AJ, Binder PS. Management of corneal ectasia after LASIK with combined, same-day, topography-guided partial transepithelial PRK and collagen cross-linking: the Athens Protocol. J Refract Surg. 2011;27(5):323-331.
4. Coskunseven E, Jankov MR, Grentzelos MA, et al. Topography-guided transepithelial PRK after intracorneal ring segments implantation and corneal collagen CXL in a three-step procedure for keratoconus. J Refract Surg. 2013;29(1):54-58.
5. Anera RG, et al. Changes in corneal asphericity after laser refractive surgery, including reflection losses and nonnormal incidence upon the anterior cornea. Opt Lett. 2003;28:417-419.
6. Cummings A. Innovations in excimer laser refractive technology – focus on the WaveLight® EX500 Excimer Laser System. Eur Ophthalmic Rev. 2010;4:44-46.
7. Data on File / WaveLight® FS200 Femtosecond Laser User Manual.
8. Stulting RD, Fant BS. Results of topography-guided laser in situ keratomileusis custom ablation treatment with a refractive excimer laser. J Cataract Refract Surg. 2016;42;11-18.
9. PERS - Declaration of Conformity for WaveLight EX500 (Nov 2017)
10. Procedure Manual EX500 (1016) rev04 2017-02-27
11. Analysis of ethanol effects on corneal epithelium PubMed: <https://www.ncbi.nlm.nih.gov/pubmed/23674759> [1]
12. Epi-LASIK: comparative histological evaluation of mechanical and alcohol-assisted epithelial separation PubMed: <https://www.ncbi.nlm.nih.gov/pubmed/?term=Epi-LASIK%3A+comparative+histo...> [2]
13. Single-Step Transepithelial PRK vs Alcohol-Assisted PRK in Myopia and Compound Myopic Astigmatism Correction PubMed: <https://www.ncbi.nlm.nih.gov/pubmed/?term=3%09Single-Step+Transepithelial+PRK+vs+Alcohol-Assisted+PRK+in+Myopia+and+Compound+Myopic+Astigmatism+Correction> [3]

---

**Source URL:** <https://www.wavelight.de/node/14511>

#### Links

[1] <https://www.ncbi.nlm.nih.gov/pubmed/23674759>

[2] <https://www.ncbi.nlm.nih.gov/pubmed/?term=Epi-LASIK%3A+comparative+histological+evaluation+of+mechanical+and+alcohol-assisted+epithelial+separation>

[3] <https://www.ncbi.nlm.nih.gov/pubmed/?term=3>